

Dear Sir/Madam,

It is nice to hear that you are interested in SeffShoe Footwear.

SeffShoe Footwear is provided after getting maximum information on the condition of your feet and Gait. Your daily walking and standing pattern and the surface on which you walk in and out of House.

Please furnish the details. Try to get the photograph as described and take the foot map of both the feet by following the instructions religiously. Try to get other medical information from your physician.

The footwear type you should wear depends upon this information and the neuropathic status of your feet. We will help you in the same. If there is any wound or callus on your feet encircle it with lipstick before taking the foot map. At our Foot Gait Laboratories we do Videography photography etc. Please note that, more the info you provide more the footwear will be of help.

Please send the form, Foot Map and photographs to us by mail, E-mail, or courier. type of footwear, the colour choice and send it to us along with the money. or deposit the amount by cheque in our account and send the slip in original to us. Keep the Xerox with you.

You can view the photographs on www.SeffShoe.com

For any help please feel free to contact us.

With Regards.

D. K. Kelkar

For Your Feet Only

Dhansai Laboratory / SeffShoe Diabetic Footwear

Address: 3,A, Uttung co-op Society, behind Bharat Bank,
off sonawala Road, Goregaon East,
Mumbai 400 063

contacts: www.dhansai.com www.seffshoe.com dhansai@gmail.com

phones: 022 26863711, 65065162 Mo. 9819526756, 9833410235

Take the Photographs like this from front and back while standing normally and looking straight ahead

If you could take video at this height while taking a few steps to and fro it will be a great help.



Slipons for Home and limited outdoor use
Rubber sole for Home use and High Grip for those whose foot is not at risk.



Sandals for Home and limited outdoor use
Rubber sole for Home use and High Grip for those whose foot is not at risk.



Diabetic Sandals for outdoor use
High Grip Sole for those whose foot **at risk**.



Shoes for outdoor use
High Grip Sole for those whose foot **at risk**.



Shoes and boots for outdoor and indoor use
High Grip Sole for those whose foot **at high risk or with Proprioception loss**.



Designer Footwear



SeffShoe Diabetic Footwear

How to take the measurements.

Seffshoes are supplied as per the size of the foot. The shoe size like 8, 9, may not match exactly with the foot. The foot measurements give much more information about one's feet. So taking precise measurement of each foot is important. follow the steps given below to take the measurements and send it along with the form and money to the following address and get your footwear as per your needs.

Material required for the measurements:-

- 1)Two sheets of plane paper, letter or legal size, depending upon the foot length.
- 2)A metal refill (Not Plastic refill or Jotter refill). This is the thin refill like plastic one. The plastic one bends and the measurement goes wrong. So use a metal refill.
- 3)If the foot problem is more serious, consult our helpline@dhansai.Com

Step 1: Put one paper on plane ground surface and keep your right foot on it.



Figure 1

Step 2: Stand normally with second foot alongside the first keeping a distance of about one foot as per figure 1

Step 3: Ask Somebody to hold the refill as shown in the figure 2. keep it vertical and draw the foot map.



Figure 2



Figure 3

Step 4: one may have to move from front to back to take the measurement. Be careful while taking the refill around heel, big toe, and near ankle bones. The fingers holding the refill should not touch the ankle bone.

Refer to figure 3 & 4



Figure 4



Figure 5

Step 5: Draw the arch as shown in figure 5 by holding the refill at about 45 degrees.

The Foot map should look similar to figure 6




Figure 6

Repeat the steps for left foot.

Write your name and address, phone No. , E-Mail etc on both the sheets.

Chose the footwear. take the photograph as described and send all the information to us by post or courier.